


NOAC



**Role of the Clinical Laboratory in Patients
receiving New Oral Anticoagulants**

NOACs

Target a single coagulation factor



Potential indications for NOACs

VTE Prophylaxis in Major Orthopedic Surgery (MOS)

Stroke Prevention in Atrial Fibrillation (SPAF)


VTE Treatment (Rx)

VTE Prophylaxis in Medical Patients (MED)

VTE recurrence prophylaxis (REC)

Acute Coronary Syndromes (ACS)

Prosthetic Heart Valves (PHV)



NOACs

What are they approved for? _____

<i>Europe</i>	MOS	SPAF	Rx	MED	REC	ACS	PHV
dabigatran etexilate Pradaxa[®]	✓	✓					X
rivaroxaban Xarelto[®]	✓	✓	✓		✓	✓	
apixaban Eliquis[®]	✓	✓				X	
edoxaban Lixiana[®]	✓ <i>Only in Japan</i>						

As of April 2013, 5th

Situations requiring laboratory control

Oral direct inhibitors of factor IIa and Xa

- **Bleeding**
- Before surgery or invasive procedure
 - the patient has taken the drug
 - in the previous 24h
 - or longer if creatinine clearance is < 50 mL / min
- Identification of sub- and supratherapeutic levels in patients
 - taking other drugs known to significantly affect pharmacokinetics
 - at extremes of body weight
 - with deteriorating renal function
- Suspicion of overdose

French GIHP* recommendation

Perioperative situations
Can I perform surgery on my patient under rivaroxaban or dabigatran?

Drug concentration

< 30 ng/mL

Perform surgery



30 - 200 ng/mL

- Wait 12h and **repeat drug measurement**
- If surgery cannot wait, **perform surgery**
- In case of **abnormal bleeding**, consider antidote (PCC or FEIBA)**

200 - 400 ng/mL

- Wait 12-24h and **repeat drug measurement**
- If surgery cannot wait, **delay surgery** as much as possible
- In case of **abnormal bleeding**, consider antidote (PCC or FEIBA)**

> 400 ng/mL

Overdose, risk of major bleeding



Pernot G, Abaladejo P.: GIHP recommendation, February 2013, available at eurekapro.fr

* Groupement d'Intérêt en Hémostase Périopératoire

** PCC: Prothrombin Complex Concentrate - FEIBA: Factor Eight Inhibitor Bypass Activity

Limitations of "non specific test"

Example of PT test

- **PT is not specific**
 - may be prolonged in case of factor deficiency, LA...
 - risk of drug concentration overestimation
- **PT baseline is not always available**
- **Apixaban induces no / limited PT prolongation**
- **PT prolongation depends on reagent sensitivity**

→ **PT is not recommended to estimate NOAC blood concentration**

Limitations of “non specific test”

Example of aPTT test

- **aPTT:** measures activity of various factors
- **Anti-IIa drug (dabigatran)**
 - curvilinear dose-response
 - steep increase at low concentrations
 - normal aPTT cannot exclude significant anticoagulant level onboard
- **Anti-Xa drugs: aPTT is not sensitive enough**
 - no prolongation or minimal prolongation...

→ **aPTT** test: “not suitable” for NOAC activity determination

Anticoagulant line

A comprehensive product portfolio

Prevention & treatment
of venous thromboembolic disease

VKA




Heparins

- UFH 
- LMWH 




fondaparinux

- calibrator 

Anti-Xa

- rivaroxaban 
- apixaban 
- edoxaban 

DTI

- argatroban 
- dabigatran 
- bivalirudin 

treatments
Stago ranges

New oral anticoagulants

PT/INR




Anti-Xa range
With calibrators and specific QC

ECA range

Validation study - Results

LCMS versus STA[®]-Liquid Anti-Xa – rivaroxaban

- Working range: 25-500 ng/mL



Slope	1,09
Intercept	-13,08
Coefficient of correlation	0,995


→ Very good correlation between LCMS and STA[®] - Liquid Anti-Xa

Validation study - Results

LCMS versus STA[®]-Liquid Anti-Xa - apixaban

- Working range: 20-800 ng/mL


Slope	0.9789
Intercept	+2.8089
Coefficient of correlation	0.994



➔ Very good correlation between LCMS and STA[®] - Liquid Anti-Xa

Anti-IIa range: ECA line

Principle



- **Chromogenic** assay derived from the Ecarin Clotting Time (ECT)
- allows the **exact DTIs concentration measurement**
- **No influence** of coagulation factors or inhibitors of the sample

Ecarin Clotting Time (ECT)

vs Ecarin Chromogenic Assay (ECA)

ECT

- Clotting test
- Reference assay, but sensitive to
 - Factor II
 - Fibrinogen level
- Not standardized
- Not commercially available

ECA

- Chromogenic assay
- NOT sensitive to factor levels
- Standardized, fully automated
- CE marked

Conclusion

- NOAC do not require routine monitoring
 - but **measurement can be helpful** in special populations and/or in special clinical situations
 - **PT and APTT are not relevant** for accurate drug concentration assessment
- Stago offers a **comprehensive portfolio** for NOAC plasma concentration measurement
 - using **specific assays** along with **dedicated calibrators and controls** (STA[®]-Liquid Anti-Xa + dedicated Cal & Ctrl / ECA line + dedicated Cal & Ctrl)